

DISTAIR S



Composition:

| | |
|--------------------------------|----------------------------------|
| PYRETHRUM EXTRACT AT 50% | 2,0 % (1,0 % of pure pyrethrins) |
| PIPERONYL BUTOXIDE..... | 10,0 % |
| COFORMULANTS UP TO..... | 100 % |

FEATURES AND USES: DISTAIR S is an insecticide pyrethrum based, to use with special electronic automatic dispensers. It is a useful remedy against flying insects, it has an immediate knock-down effect of flies, mosquitoes, tiger mosquitoes, moths, etc. DISTAIR S is suitable for the disinfection of domestic and civil (eg. Housing, living rooms, toilets, meeting rooms, schools, communities, canteens, shops, workshops, warehouses, food industries, hotels, restaurants, bars, barracks, stables, animal shelters, transportation such as ships, aircraft, trains, buses, motor homes, etc.). The repellent effect of pyrethrum extract can help and prevent the infestation of the premises treated by insects.

HOW TO USE: Insert the can into the dispenser and activate it following the instructions in the package. For best results put the electronic dispenser securing it to a wall or resting it on a mobile at least two meters from the floor and 30 cm from the ceiling. DISTAIR S is released at programmable intervals from 7.5 to 30 minutes in very small quantities, ensuring a constant protection from insects. A can of 250 ml lasts up to 4 weeks with a use for 24 hours per day at an interval time of 15 minutes. The system is designed for the control of insects in environments with dimensions up to 30 m². DISTAIR S can also be used manually: keep the can upright, press the nozzle and spray upwards towards the center of the room, moving in short spurts. For a medium sized environment does not exceed a vaporization of 5-6 seconds.

PRECAUTIONS FOR USE: Do not use in places where there are aquariums. The non-covered foods, prepared or in preparation, should be located more than 3 meters away from the can.

WARNING: For proper use of this product, follow the instructions given on the label .

NOTES: CONFIDENTIAL INFORMATION TO DOCTORS, HYGIENISTS, VETERINARY AND PARASITOLOGISTS.
Volume 250 ml - Net weight g . 160-3000 splashes - dispensing valve dosage from 75 microliters

Code
028



ORMA S.r.l. - Via U. Saba, 4 - 10028 Trofarello (To) Italy
TEL. +39 011.64.99.064 - FAX +39 011.68.04.10
www.ormatorino.com - e-mail: aircontrol@ormatorino.it